



# Chronic Maintenance In -Center Hemodialysis Standing Orders (Addendum)

Chronic Maintenance In -Center Hemodialysis Standing Orders ± Medications PRN ,  
item 11b clarification:

11b.

### Staff Application

Apply a moderate layer of lidocaine 2.5% / prilocaine 2.5% cream (or other generic equivalent) at the site of needle placement 15-60 minutes (preferably 30 or greater minutes) prior to start of hemodialysis.

### Patient Application

Dispense to patient: One 30-g ram tube of lidocaine 2.5% /prilocaine 2.5% .

Directions for use: Apply a moderate layer (approximately 2.5 grams) of cream at site of needle placement 15-60 minutes (preferably 30 minutes or greater) prior to start of treatments. Not to exceed one 30-g ram tube per 12 dialysis treatments.

ODWWKHZ MDYDUD  
Physician Name (Please Print)

\_\_\_\_\_

\_\_\_\_\_  
Physician signature  
(see Initial Orders)

\_\_\_\_\_  
Date

Patient Name : \_\_\_\_\_ NKC# \_\_\_\_\_